

FABT[®]

PRODUCT CATALOGUE

DEFIBRILLATION
EURODEFIPADS[®]



†EURODEFIPADS®†

MULTIFUNCTION PADS



High quality materials, continuous process control and extreme care in assessing the compatibility and the conformity to the international standards make EURODEFIPADS® excellent and safe devices. Easy to be used, they assure optimal clinical performances both in emergency and elective cases.

EURODEFIPADS® is a complete line of adult and pediatric disposable multifunction electrodes for defibrillation, synchronized cardioversion, external cardiac stimulation and ECG monitoring with different connectors compatible with the main defibrillators present in the market.

Radiotransparent models (with radiolucent cable and pads) have been designed for use in EP labs, haemodynamics and during radiographic and angiographic procedures. Lead-out preconnectable models (marked with suffix W) have cable and terminal connector outside the pouch, allowing to connect pads to the defibrillator before opening pouch; they have been intended especially for emergency situations.

Euro Defi Pads are compatible with monophasic and biphasic defibrillators (manual, semi-automatic and automatic) and they comply with **European Directive and International Standards** (IEC/EN 60601-1, IEC/EN 60601-2-4, ANSI/AAMI DF-80, EN ISO 10993-1)

ADULT



PAEDIATRIC



RADIOTRANSSPARENT
























LEAD OUT








Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model. Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

		ADULT		PAEDIATRIC				
		Radiopaque	Radiotransparent	Radiopaque	AED	Radiotransparent		
<p>CONNECTORS</p>   		<p>BRANDS</p> <p>Philips Medical® Laerdal Medical® Agilent Philips Medical</p>		 F7950 F7950W F7950U F7950UW	 F7750 F7750W	 F7950P F7950PW F7950U F7950UW	 F7950PL F7950WPL	 F7750P
		<p>Zoll Medical Corp® M&B</p>		 F7951 F7951W	 F7751 F7751W	 F7951P F7951PW	 	 F7751P
<p>F7952PL F7952WPL</p>   <p>other</p>		<p>Physiocontrol® Bexen Cardio Cardioline® Mindray Comen</p>		 F7952 F7952W F7952/200	 F7752 F7752W	 F7952P F7952PW	 F7952PL F7952WPL	 F7752P F7752PW

EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model. Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

		ADULT		PAEDIATRIC		
		Radiopaque	Radiotransparent	Radiopaque	AED	Radiotrasparent
CONNECTORS	BRANDS					
	Agilent Laerdal Medical® Hewlett-Packard	F7953 F7953W	F7753	F7953P F7953PW		F7753P
	Schiller Esaote	F7954 F7954W	F7754	F7954P		F7754P
	Nihon Kohden Welch Allyn® Innomed MRL	F7955 F7955W	F7755 F7755W	F7955P F7955PW		F7755P
	Schiller	F7956 F7956W	F7756	F7956P F7956PW		F7756P
	Cardiac Science	F7957W			F7957PL	
	CU Medical Systems, inc	F7958 F7958W		F7958P F7958PW	F7958PL	
	CU Medical Systems, inc Cmos Drake	F7959 F7959W	F7759	F7959P	F7959PL F7959WPL	
	Metrax Primedic	F7960W	F7760W			
	Nihon Kohden	F7961U F7961UW F7961UWXL	F7761	F7961U F7961UW F7961UWXL F7961P F7961PW		F7761P
	Marquette GE® Esaote	F7962	F7762	F7962P		
	Nihon Kohden	F7963				
	Weinmann Cardiaid	F7964W		F7964PW		

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model. Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

		ADULT		PAEDIATRIC		
		Radiopaque	Radiotransparent	Radiopaque	AED	Radiotrasparent
CONNECTORS	BRANDS					
	Corpuls	F7965 F7965/PLUS		F7965P		
F7966W  F7966WPL 	Defibtech	F7966W			F7966WPL	
F7967  F7967W  F7967P  F7967PW 	Schiller	F7967 F7967W		F7967P F7967PW		
F7969W  F7969PW 	Mediana	F7969W		F7969PW		
 	GE®	F7981		F7981P		
	Instramed Aed	F7984W		F7984PW		
	Mediana Tecnogaz	F7986 F7986W F7986UW		F7986P F7986PW F7986UW		
	Zoll Medical Corp®	F7987W		F7987PW		
	Nihon Kohden	F7989UW		F7989UW		
Not for sale in France, Germany and UK						
	Laerdal Medical®	F7970		F7970P		
	Physiocontrol®	F7971				
	Marquette GE® Esaote	F7972		F7972P		

EURODEFIPADS®




STIMULATION PADS

CE 2797

Besides EURODEFIPADS® models, some additional pads have been developed only for external cardiac stimulation.



F7955/STIM

CONNECTORS	BRANDS		
	Zoll Medical Corp® (M series, pacer version only)	adult paediatric	F7951/STIM F7951P/STIM
	Physiocontrol®	adult	F7973
	Nihon Kohden Innomed	adult	F7955/STIM

TRAINER PADS

CE

Defibrillation training pads (adult and paediatric) are equipped with snap connection and must be used with the appropriate connection cable. The pad's special shape makes the application extremely easy both on manikins and on other surfaces.

PADS

F7999	adult
F7999P	paediatric



F7999P



F7999



F5352

CONNECTING CABLES

CONNECTORS	BRANDS	
	Philips Medical® Laerdal Medical® Agilent Philips Medical	F5350
	Physiocontrol® Osatu Bexen Cardio Cardioline® - Mindray	F5352
	Welch Allyn	F5355
	Schiller	F5356
	Weinmann Cardiaid	F5364
	Defibtech	F5366
	Mediana	F5369
	Mediana Tecnogaz	F5386
	Nihon Kohden	F5389



FIAB SpA
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy
Tel. +39 0558497999 - Fax +39 055844381
export@fiab.it - www.fiab.it

